

## HEINE UniSpec® Disposable Ear Specula



2,5 mm dia.



4 mm dia.



- ∴ **Disposable.** Hygienic and safe.
- ∴ **Original HEINE design and manufacture.** Fit positively to the otoscope.
- ∴ **Assured quality.** No sharp edges.
- ∴ **A clean, new speculum for each and every patient.** The right size visible and instantly available.

### UniSpec Disposable Ear Specula (grey)

#### for BETA100, K100, Operating Otoscope

Pack of 1000, 2,5 mm dia. (infants) 1000 pcs./1 pack	<b>B-000.11.242</b>
Pack of 1000, 4 mm dia. (adults) 1000 pcs./1 pack	<b>B-000.11.241</b>
UniSpec tip dispenser	<b>B-000.11.146</b>

⊗: Single use only.

## Reusable Specula



2.2 mm dia.



2.8 mm dia.



3.5 mm dia.



4.5 mm dia.



5.5 mm dia.



10 mm dia.

- ∴ **Moulded in high-density SANALON S, smooth surface.** Easy to clean, can be disinfected with conventional solutions. Autoclavable at 134 °C.
- ∴ **Ergonomic Design.** Painless and atraumatic insertion.
- ∴ **Metal security lock with bayonet fitting.** Long-lasting, positive fit to the otoscope.

### Reusable Specula (black)

#### for BETA100, K100 and Operating Otoscope

1 set = 5 ear specula and 1 nose speculum	<b>B-000.11.221</b>
3 specula 3.5/4.5/5.5 mm dia.	<b>B-000.11.222</b>
Ear speculum 2.2 mm dia.	<b>B-000.11.215</b>
Ear speculum 2.8 mm dia.	<b>B-000.11.216</b>
Ear speculum 3.5 mm dia.	<b>B-000.11.217</b>
Ear speculum 4.5 mm dia.	<b>B-000.11.218</b>
Ear speculum 5.5 mm dia.	<b>B-000.11.219</b>
Nose speculum 10 mm dia.	<b>B-000.11.220</b>