

Reusable VET specula

for BETA 200 VET or BETA 100 VET with adaptor



- : Made of high-density SANALONS with a smooth external surface. Dust and germs are easily removed. Disinfection and autoclaving at up to 134 °C is possible.
- : Ergonomic shape. Atraumatic.
- : Convenient fixation system. Positive fit to otoscope head.

Reusable SANALON S specula	
1 set = 3 reusable SANALON S specula	G-000.21.213
61 mm long / 4 mm dia.	G-000.21.210
61 mm long / 6mm dia.	G-000.21.211
61 mm long / 9 mm dia.	G-000.21.212

HEINE AllSpec® Disposable tips

for BETA100 VET with adaptor



- : For single use. Hygienic, no cross-infection.
- : Original HEINE design. Firm fit to the otoscope head.
- : Assured quality. No sharp edges.

Ch.	A	F
11/2	ASSONE"	C
[01]		





AllSpec Disposable tips (grey)	
Pack of 1000 pcs. 2.5 mm dia. [03]	B-000.11.128
Pack of 1000 pcs. 4mm dia. [03]	B-000.11.127
Case of 10 dispenser packs of 250 tips, 2.5 mm dia [01]	B-000.11.151
Case of 10 dispenser packs of 250 tips, 4mm dia [01]	B-000.11.150
AllSpec Tip-dispenser [02]	B-000.11.148

AllSpec Tip-dispenser [02]

The right size, visible and instantly available. A clean new tip for each and every patient. Contents: 50 pcs. Ø 2.5 mm and 50 pcs. Ø 4 mm.

(X): Single use only.

Insufflation Bulb



For BETA 200 VET and BETA 100 VET Otoscope.

: For pneumatic testing.

Insufflation bulb	B-000.11.240
-------------------	--------------